



Name: _____

Class: _____

Experiment worksheet

4.2 Acid reactions depend on strength and concentration

Pages 92–93 and 203

Experiment 4.2: Acid titrations

Aim

To compare the reactions of a strong acid, hydrochloric acid, and a weak acid, ethanoic acid, (common name acetic acid).

Materials

- Dropper bottles containing:
 - 0.1 M hydrochloric acid (HCl)
 - 0.1 M ethanoic acid (acetic acid) (CH_3COOH)
 - 0.1 M sodium hydroxide (NaOH)
 - 1 M hydrochloric acid (HCl)
 - 1 M ethanoic acid (acetic acid) (CH_3COOH)
 - universal indicator solution
- pH colour chart
- Small pieces of magnesium ribbon
- 4 test tubes and test-tube rack
- Dropping pipette
- Matches



CAUTION: Wear a laboratory coat, gloves and safety glasses throughout this experiment. Most of the substances used in this experiment need to be handled with care.

Class:

PART A

- 1 Draw up a table to record each test and the results for each acid.
- 2 Place 2 mL of 0.1 M hydrochloric acid in one test tube and add 2 drops of universal indicator solution. Record the colour of the indicator and the corresponding pH from the colour chart.
- 3 Repeat step 2 with 0.1 M ethanoic acid, using a fresh test tube.
- 4 To the first test tube add 0.1 M sodium hydroxide drop by drop, counting the drops, until the solution is neutral (i.e. pH = 7).
- 5 Repeat step 4 with ethanoic acid.

PART B

- 1 Add 2 mL of 1 M hydrochloric acid to a fresh test tube.
- 2 Add a small piece of magnesium ribbon to the test tube and invert a clean test tube over the top so there is only a small gap between the test tubes.
- 3 Record your observations.
- 4 Lightly touch the base of the bottom test tube. Record your observations of the temperature of the mixture.
- 5 When the reaction has ceased, light a match and hold it just inside the inverted test tube. Do you hear a loud popping sound? This is evidence of hydrogen gas being produced.
- 6 Repeat steps 1–6 with 2 mL of 1 M ethanoic acid.

Results

Record your results in an appropriate table.

[illegible]



Name: _____

Class: _____

Discussion

- 1 When you tested the pH of the two acids, you used the same concentration (0.1 M).
 - a Why were they compared at the same concentration?

 - b Why did they have a different pH?

 - c What can be concluded about the strength of ethanoic acid compared with the strength of hydrochloric acid?

- 2 Compare the number of drops of sodium hydroxide used to neutralise each acid. Is this what you expected? Explain using your results.

- 3 Write a balanced equation for each neutralisation reaction.

- 4 The pop test is the standard test for hydrogen gas. The 'pop' sound is a mini-explosion due to the combustion of hydrogen gas in air, which is a very exothermic (heat producing) reaction.

The equation for the reaction is: $2\text{H}_2(\text{g}) + \text{O}_2(\text{g}) \rightarrow 2\text{H}_2\text{O}(\text{l}) + \text{energy}$

- a Did your tests confirm that hydrogen gas was produced?

 - b Was there a difference in how fast the reactions with the two different acids occurred? If so, suggest why.

- 5 Write a balance chemical equation for the reaction between the two acids and the magnesium ribbon.



Name: _____

Class: _____

Conclusion

1 What do you know about neutralisation reactions?

2 What do you know about reactions between metals and acids?

3 What do you know about the difference between strength and concentration of acids?
