# Le Châtelier's Principle

#### INTRODUCTION

Le Châtelier's Principle states that when stress is placed on a system in equilibrium, the system will react to relieve the stress.

*Stress* refers to changes in concentration, pressure, volume, and temperature. Consider the exothermic chemical reaction:

$$SO_2(g) + 0.5 O_2(g) = SO_3(g)$$

**Concentration:** If the reaction has reached equilibrium, and we disturb this reaction at equilibrium by adding  $SO_2$ , the added  $SO_2$  would be reduced in quantity by the reaction moving to the right, forming more of the  $SO_3$ .

**Pressure:** For this reaction, an increase in pressure on the system would also move the reaction to the right. There are 1.5 total moles of gas on the left of the arrows, and only one mole of gas on the right. As the reaction moves to the right, the number of moles present would reduce, and so the pressure would reduce.

**Volume:** For this reaction, if the volume of the container is increased, the reaction would move to the left, where there are more moles of gas than there are on the right. Increases of volume and decreases of pressure are complementary processes. They can be effected by a piston-cylinder reaction vessel.

**Temperature:** The reaction is exothermic. It gives off heat. If the temperature on the system were increased, the reaction would move to the left, thereby absorbing some of the heat added.

To easily apply Le Châtelier's principle, think of the arrows as a funnel for the reaction. The material on the left side can flow to the right, and the right side can flow to the left. If you "push" or "pull" on one side of the reaction by changing the concentration of a reactant, the quantities will "flow" through this funnel in a way to offset the change. If you increase the pressure (or decrease the volume) for the system, the reaction will flow in a direction to offset the increase, if possible. If you change the temperature, the system will flow in the direction that offsets the temperature change.

#### **EXPERIMENT**

## **Supplies**

You will use test tubes from your drawer. You will estimate quantities of reactants, adding solutions from transfer pipets attached to the reagent bottles. These are marked to the proper volume. You will add other reactants by drop.

Except for Part 6, all of the reactions shown take place in water. The (aq) appendage, which shows that a substance is in water, is left off so that the equations would fit in the space allowed.

#### Part 1: Chromate/Dichromate Equilibrium.

$$2\text{CrO}_4^{2-} + 2\text{H}_3\text{O}^+ \rightleftharpoons \text{Cr}_2\text{O}_7^{2-} + 3\text{H}_2\text{O}$$
yellow orange

Place about 1 milliliter of 1 M K<sub>2</sub>CrO<sub>4</sub> in a test tube. Add several drops of 3 M H<sub>2</sub>SO<sub>4</sub>. Stir. Now add drops of 6 M NaOH with stirring until a change takes place. Add more H<sub>2</sub>SO<sub>4</sub>. Record your observations and explain what happened using Le Châtelier's principle. You need to know that H<sub>2</sub>SO<sub>4</sub> produces H<sub>3</sub>O<sup>+</sup>ions in solution, and NaOH produces OH<sup>-</sup>ions which remove H<sub>3</sub>O<sup>+</sup>ions.

Dispose of the solution in the waste container marked "Chromate Waste".

### Part 2: Indicator Equilibria.

There is a class of dyes that changes color depending on the concentration of H<sub>3</sub>O<sup>+</sup> ions in solution. These dyes are called acid-base indicators. The generic reaction

$$HIn + H_2O \rightleftharpoons H_3O^+ + In^-$$
one color another color

shows how the indicators work. They are usually carbon containing compounds with benzene rings. The concentration of H<sub>3</sub>O<sup>+</sup> required to shift the equilibrium differs for different indicators. Methyl orange is red in the HIn state, yellow-orange in the In<sup>-</sup> state. If the H<sub>3</sub>O<sup>+</sup> is greater than 10<sup>-4</sup> (pH of 4), methyl orange shifts to the HIn state. Phenolphthalein is colorless in the HIn state, pink in the In<sup>-</sup> state. If the H<sub>3</sub>O<sup>+</sup> is greater than 10<sup>-9</sup> (pH of 9), phenolphthalein shifts to the HIn state.

Place about 1 milliliter of deionized water into each of two test tubes. Add a drop of methyl orange to one, a drop of

1

phenolphthalein to the other. Add one drop of 6 M HCl to each test tube. Add two drops of 6 M NaOH to each test tube. Record your observations and explain what happened using Le Châtelier's principle. Here is the generic indicator equation rewritten for methyl orange:

$$HMO + H_2O \rightleftharpoons H_3O^+ + MO^-$$
  
red yellow-orange

HCl increases the H<sub>3</sub>O<sup>+</sup>, and NaOH decreases the H<sub>3</sub>O<sup>+</sup>. Rewrite the generic equation in a similar way for phenolphthalein in the explanation section on the next page.

Dispose of the solutions in the sink.

#### Part 3: Acetic Acid/Acetate Equilibrium.

$$HC_2H_3O_2 + H_2O = H_3O^+ + C_2H_3O_2^-$$

Pour about 1 milliliter of 0.1 M  $HC_2H_3O_2$  into each of two test tubes. Add a drop of methyl orange solution to each test tube. Add 1 M  $NaC_2H_3O_2$  solution to one test tube a drop at a time until there is a difference in color between the two solutions. Record your observations and explain what happened using Le Châtelier's principle. When sodium acetate,  $NaC_2H_3O_2$ , dissolves in water, it breaks up into the  $Na^+$  and the  $C_2H_3O_2^-$  ions. The methyl orange changes color with changing  $H_3O^+$  concentration, so refer to Part 2 of the experiment to help you interpret what the addition of acetate ions did to the equilibrium.

Dispose of the solutions in the sink.

## Part 4: Ammonia/Ammonium Equilibrium.

$$NH_3 + H_2O = NH_4^+ + OH^-$$

The OH $^{-}$  ion concentration controls the  $\mathrm{H_3O^+}$  ion concentration in water solutions. As the OH $^{-}$  ion concentration gets larger, the  $\mathrm{H_3O^+}$  ion concentration gets smaller, and vice versa. The [OH $^{-}$ ] in the original aqueous solution of ammonia, NH $_3$ , is about  $10^{-3}$ , which makes the [H $_3\mathrm{O^+}$ ] about  $10^{-11}$ .

Pour about 1 milliliter of 0.1 M NH<sub>3</sub> into each of two test tubes. Add a drop of phenolphthalein solution to each test tube. Add 1 M ammonium chloride, NH<sub>4</sub>Cl, solution to one test tube a drop at a time until there is a difference in color between the two solutions. Record your observations and explain what happened using Le Châtelier's principle. When NH<sub>4</sub>Cl dissolves in water, it breaks up into the NH<sub>4</sub><sup>+</sup> and the Cl<sup>-</sup> ions. As you were instructed in Part 3, use the information from Part 2 to help you interpret what the addition of ammonium ions did to the equilibrium.

Dispose of the solutions in the sink.

#### Part 5: Iron/Thiocyanate Equilibrium.

$$Fe^{3+} + SCN^{-} \rightleftharpoons Fe(SCN)^{2+}$$
pale yellow deep red

Pour about 1 milliliter of deionized water to each of three test tubes. Add 1 drop of 0.1 M Fe( $NO_3$ )<sub>3</sub> and 1 drop of 0.1 M KSCN to each test tube. Add two more drops of 0.1 M Fe( $NO_3$ )<sub>3</sub> to one of the tubes, and two more drops of 0.1 M KSCN to another of the test tubes. Try to compare the colors of the solutions as quantitatively as possible. Record your observations and explain what happened using Le Châtelier's principle.

Dispose of the solutions in the sink.

# Part 6: Temperature Effect.

$$CoCl_4^{2-} + 6 H_2O \rightleftharpoons Co(H_2O)_6^{2+} + 4 Cl^- + energy$$
  
purple pink

Pour about 2 milliliters of 0.15 M CoCl<sub>2</sub> (in methanol) into a test tube. Add deionized water drop by drop with stirring until the blue color *just changes* to pink. Do not add excessive water. Place the test tube in a water bath in the fume hood. The temperature in the water bath should be 65 to 70 °C. When the color changes, remove the test tube from the hot water bath, and cool it by swirling it in tap water. (If you do not see a color change, you may have added too much water. Remake the solution. Take care to follow directions.) Record your observations and explain what happened using Le Châtelier's principle.

Dispose of the solution in the waste container marked "Cobalt Waste".

#### Part 7: Saturated Equilibrium.

$$NaCl(s) = Na^{+}(aq) + Cl^{-}(aq)$$

Pour about 1 milliliter of saturated NaCl solution into a test tube. Add 12 M HCl drop by drop with stirring until you see a change taking place in the test tube. Continue to add 3 or 4 more drops, pausing after each drop. Record your observations and explain what happened using Le Châtelier's principle.

Dispose of the solution in the sink.

OBSERVATION TABLE		
R ewrite the chemical equation for each part in the explanation .		
Part 1 Observations		
Explanation		
Part 2 Observations		
Explanation		
Part 3 Observations		
Explanation		
Part 4 Observations		
Explanation		
Part 5 Observations		
Explanation		
Part 6 Observations		
Explanation		
Part 7 Observations		
Explanation		

Name	Grade	Date

## **QUESTIONS**

1.	What you did in Parts 3, 4, 5, and 7 are exam—les of the <i>common ion effect</i> . Why do you suppose it has that name?
	Give an explanation of the common ion effect based on what you did and what you observed.

2. Rewrite the chemical equation from Part 6 using 
$$\Delta H$$
 notation for energy. Simply indicate if the  $\Delta H$  is + or -.

3. Write the following equation in the manner the equation in Part 6 is written (note the energy term): 
$$2 \text{ H}_2\text{O}(1) + \text{O}_2(g) \quad 2 \text{ H}_2\text{O}_2(aq) \qquad \Delta H = 196.1 \text{ kJ}$$

