10 GENERAL SCIENCE BIOLOGY INVESTIGATION: DNA EXTRACTION 2016

Name:		Teacher: Miss Cerny
Form:		Due date:
	1/2.1	

Living things are made up of cells which contain the genetic code which distinguishes them from other living thing. This code is found in the chemical inside the nucleus of cells and is known as DNA.

Plagiarism

You must write in your <u>own words</u>, not copy sentences word for word from another student or another source.

MARKING VEY

Plagiarising = instant zero on assignment.

Assessment policy

Extensions are not given on the day an assessment is due.

Give me a sick note/legitimate reason from parent BEFORE due date = new negotiated due date.

One day late = -20% taken off mark.

Two days late = -40% taken off mark.

Three days late = mark of zero given and students are required to attend a detention and are still required to submit the assignment.

If you are not at school the day this assignment is due, please email it to me by the 4pm deadline.

s.cerny@aranmore.wa.edu.au

MARKINGKET

Aim: to extract DNA from a strawberry.

Introduction

Write an introduction using $\underline{\textbf{full sentence paragraphs}}$ that includes the following points.

◆ Describe three reasons why scientists would need to extract DNA from living things today. (3	marks)
Friedrich Miescher (1)	
In 1869	
Germany ()	
-Extract DNA to study genes involved in cancer.	
- Extracted DNA can test for genetic diseases	
in neuborn babies.	
- To malyse forensic material in criminal	
investigations	
- To accurately organise/sort organisms into classes	5
-To allow for the genetic modification of	
organisms.	
Tope mark for each three	
Take a mark off if not	
written in Jul sentences.	
	,
•	

Potential hazards

You will be using isopropyl alcohol in this investigation. Describe three potential hazards of using this chemical in the science laboratory. (3 marks)
- Hazardous in case of eye contact it irritates and can damage eyes).
- Hazardous in case of ingestion (can cause nausea, vomiting, etc.)
- Hazardous in case of inhabation (concause headaches, dizziness etc.)
- Slightly hazardors in case of skin Gntactemay initate the skin.
- Hazardous if in Ontact with open flames /sparks as it is
highly flanmable.
Anythree- describe NOT state)
- Needs to relate to science lab
Take a mark off if not
enough into (description is
Explain two things you will do to help prevent hazards occurring while using isopropyl alcohol in this investigation. Be sure to explain why you will take these precautions. (4 marks)
-wear safety glasses so isopropyl does not enter
theeyes.
- Do not crouch in Front of the bench as if
things spill it may go on you.
- Wear gloves to avoid isopropyl coming into (OR)
contact with skin
- Ventilate the rooms so that fumes are able
to dissipate.
One mark for way to prevent hazard & one mark for explanation as to why
Title: (write a title for the investigation, needs to be detailed). (1 mark)
eg. Extraction of DNA from a strawberry
eg. Isolation of Strawberry DNA to measure mass

Materials: (list all materials used, be specific with number of	fitems used and amounts/sizes). (2 marks)
100ml distilled water	1 x pair + Welzers
14 tsp salt	1x stirring rod
1 x 250 ml beaker	1x electronic balance
10ml dish soap	1x watch glass
1 x plastic zip-loch bag	1 x large measuring cylinder
1x strawberry	1x small measuring cylinder
1 x sieve	
Sml chilled isopropyl	(option, some groups may have used a plastic spoon).
Method	
1. Measure 100 ml of distilled water into a 250 ml beaker.	
2. Add ¼ tsp of salt to the water in the beaker and mix using	the stirring rod.
3. Add 10 ml of dish soap to the water and salt solution and	gently mix with the stirring rod.
4. Place one strawberry into a plastic zip-lock bag.	
5. Pour the soap, water and salt solution into the zip-lock ba	g.
6. Remove as much air from the bag as possible and seal it.	
7. Use your hands to squash and mush the strawberry inside	e the bag until there are no large pieces remaining.
8. Pour the liquid from the bag through a sieve sitting over t	he original beaker (that has been rinsed).
9. Use the stirring rod to press the strained strawberry throu	ugh the sieve.
10. Add 5 ml of chilled isopropyl alcohol to the solution (do	
11. Gently use the tweezers to remove the DNA strands from	n the beaker and observe.
Results	
Describe the DNA that you extracted from the strawberry (v	vhat does it look like). (1 mark)
eg: The DNA strands looke	d like colourless,
clear strings of salivo	d like colourless, a with bubbles attached

Work out how much of the strawberry's mass is made up of DNA. Show this a working out. (You can research how to do this). $\underline{\ell}g$	ns a percentage and show all of your (4 marks)
Mass of strauberry is: 12.85g	
Mass of strauberry is: 12.85g Mass of DNA is: 0.49g	
Amount of Stranberry's mass made up of ONA	
	mass of strauberry
Amount of strawberry's mass made up of DNA =	
	12.85
	3.87. (1)
List the steps you took to complete the above task. (The mass of the strawbe balance etc.)	erry was measured using an electronic (3 marks)
1. Strawberry was placed on electron	ic balance
and the mass was recorded.	
2. The DNA was placed in a petri	dish.
3. An empty petri dish was place	
electronic balance and re-zero	o'd
4. The petri dish with DNA was	
electronic balance & the mass was	recorded.
5. A formula was used to calculate	how much of
the strawberry's mass was made	V
OR something similar	

Discussion


State the function of the dish soap in the extraction of DNA from the strawberry.	
The dish soap helps to dissolve the cell	
membranes.	
S CONTRACTOR OF THE STATE OF TH	(1 moule)
State the function of the salt in the extraction of DNA from the strawberry.	(1 mark)
The salt breaks up protein chains that	
hold nucleic acids together.	
Explain why the chosen fruit for DNA extraction was strawberries and not a different type of fruit.	(1 mark)
The strawberry has more DNA than any	
other fruit.	

Reference list

Minimum of two references.		(1 mark)	
Referenced using	the APA referencin	g style.	(1 mark)
	If there	is even	
	01	e mistake (missing full stop, etc)	
		e mistake (missing full stop, etc) then they do not get	
		this mark.	
4-11-11-11-11-11-11-11-11-11-11-11-11-11			
Presented neatly a	and clearly.		(1 mark)
Correct grammar	and spelling.		(1 mark)
Uses scientific lan	guage.		(1 mark)
Total mark:	/32		
Percentage:	%		
rercentage.	70		
			
Teacher's comme	nts:		

APA referencing guide

Book



Book with two authors

Shermer, M., & Benjamin, A. (2006). Secrets of Mental Math: The Mathemagician's Guide to Lightning Calculation and Amazing Mental Math Tricks. New York, USA: Three Rivers Press.

Book with three or more authors

Bulliet, R. W., Alley, R. B., Broecker, W. S., & Denton, G. H. (2011). *The fate of Greenland: Lessons from abrupt climate change*. Cambridge, MA: MIT Press.

Book with editor

Kasdorf, W. E. (Ed). (2003). *The Columbia guide to digital publishing*. New York, USA: Columbia University Press.

Website




Website with no author (put title of page at the front)

Aranmore Catholic College. (2014). Retrieved from http://www.aranmore.wa.edu.au/

Website with no date

The Zadkine Museum. (n.d). Retrieved http://www.france.fr/en/museums/zadkine-museum.html

Video file



Encyclopaedia or dictionary

